

User Manual

Version: Ver 1.3

Revision date: 05-2022

Product name: Fingertip Pulse Oximeter

Doc.:IFU-FPO

Model:

YK-80A, YK-80B, YK-80C, YK-81A, YK-81B, YK-81C, YK-82A, YK-82B, YK-82C, YK-83A, YK-83B, YK-83C, YK-84A, YK-84B, YK-84C, BPX800US, YK-81CEU, BPX800CA.

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive features for artery SPO₂ and PR detection. Being portable, it is able to measure SPO₂ and PR values quickly and precisely.

General Description

Haemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO₂) that compounded with oxygen and that of all combinativable haemoglobin (Hb) in blood. In other words, it is the saturation of Oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce haemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human haemoglobin saturation. As a result, such symptoms as megrim, vomiting and asthenia might appear to patients. Hence, it is very important to know hemoglobin saturation of patient timely in clinical medical aspects.

The fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

Intended use

The fingertip pulse oximeter can be used to measure human haemoglobin saturation and pulse rate through finger , it can be used in hospitals, families, schools and medical centers.


User Group

-People who need blood oxygen measurement .

Contraindication: not found

⚠ Note ⚠:

- 1.The image in the instruction may have slight differences with the actual instruments.
- 2.Technical parameters and appearance change, without prior notice.




Product include: main machine and SPO2 sensor。

Measurement principle

The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin and oxyhemoglobin (HbO₂) in red light and infrared light zones. Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm red light and 940nm infrared light) can be focused onto human nail through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on two groups of diodes through process in electronic circuits and microprocessor.


Diagram of Operation Principle



1. Red-ray and Infrared-ray receiving diode
2. Red-ray and Infrared-ray transmitting diode

Precautions for use

1. Do not use the Fingertip Pulse Oximeter together with MRI or CT equipment.
2. Explosion hazard: Do not use the Fingertip Pulse Oximeter in an explosive atmosphere.
3. The Fingertip Pulse Oximeter is intended only as an adjunct in patient assessment. Doctors should make diagnosis in conjunction with clinical manifestation and symptoms.
4. Check the Fingertip Pulse Oximeter sensor application site frequently to make sure that the circulation and skin integrity of patient are under good condition.
5. The sensor of the oximeter is not suitable for contacting the adhesive tape, which may lead to the error of measurement data or mistaking that there are blisters on the tested skin.
6. Please read the manual carefully before your operation.
7. The Fingertip Pulse Oximeter has no SpO₂ alarm, it is not for continuous monitoring.
8. When used for a long time, it will produce discomfort or tenderness, especially for patients with microcirculation disorders. Please change the test sites periodically according to different situations of patients. **The test site must be changed and the skin integrity and circulation condition of the patient must be checked at least every 2 hours, and the correct adjustment must be made.**
9. Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.
10. Significant levels of dysfunctional hemoglobins (such as carboxyl-hemoglobin or methemoglobin) may cause inaccurate reading.
11. Intravascular dyes such as indocyanine green or methylene blue may cause inaccurate reading.
12. SpO₂ measurements may be adversely affected in the presence of high ambient light. Please




- shield the sensor area (with a surgical towel or direct sunlight, for example) if it is necessary.
13. Unexpected action may cause inaccurate reading.
 14. Medical signal with high frequency or interference caused by defibrillator may lead to inaccurate reading.
 15. Venous pulsations may cause inaccurate reading.
 16. It may cause inaccurate reading when the positions of sensor and blood pressure cuff are on the same arterial catheter or intravascular line.
 17. Hypotension, severe vasoconstriction, severe anemia, or hypothermia may cause inaccurate reading.
 18. It may cause inaccurate reading by giving use of cardiogenic to patient after his cardiac arrest or when he is in quiver.
 19. Bright nail or painted nail may cause inaccurate SpO2 reading.
 20. Do not use this product if you are allergic to ABS, black silicone pad and other materials.
 21. If the performance is inconsistent with the description or changes, stop using immediately and contact the manufacturer.
 22. Measuring function should not be used to evaluate oximeter accuracy.
 23. The effects of degraded sensors and electrodes, or loosened electrodes, that can degrade performance or cause other problems.
 24. Please place the oximeter in a place where children, pets, etc. can not be touched, so as to prevent falling, biting and affecting the product performance.
 25. The height of the oximeter should not exceed two meters.
 26. Do not use beyond the service life of the product, otherwise the accuracy of the instrument will be affected
 27. No modification of this equipment is allowed.
 28. Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm(12 inches) to any part of the oximeter, including cables specified by the manufacturer. Otherwise, degradation of the performance of this Fingertip Pulse Oximeter could result. Otherwise, degradation of the performance of this oximeter could result.
 29. Oximeter should be avoided in places with poor ventilation and high dust and lint content.
 30. Avoid direct sunlight and strong light sources.
 31. Oximeter cannot be serviced and maintained during use.
 32. The internal structure of the Fingertip Pulse Oximeter adds magnetic rings, magnetic beads, and conductive cloth to spray conductive paint to avoid electromagnetic interference, so as to prevent adverse events to patients and operators due to electromagnetic interference.
 33. The Fingertip Pulse Oximeter can be maintained and calibrated once every two years, and the basic safety and basic performance of the Fingertip Pulse Oximeter have been guaranteed
 34. It is necessary to change the placement position of the oximeter regularly, check the placement position of the oximeter every 2 hours, check the impact of the oximeter on the skin, and move the oximeter to different positions.

Features

◆ OLED display

1. Product adopts double color OLED display, can show the six different display mode.

【 If the hand movements, under the effect of gravity induction, the interface can have four




different kinds of display mode (suitable for matching gravity induction function instrument)】

2. Low-power consumption, continuously work for more than six hours with two AAA batteries
3. Low voltage indicator
4. In the absence of signals, the product will be in after 8 seconds to enter standby state.
5. Small in volume, light in weight, and convenient to carry

【 Instrument has 5s automatic signal detection function, when you insert finger, timely will automatically start;Automatic startup function instrument (applicable to Automatic startup function instrument)】

Operation Instructions

1. Install two AAA batteries into battery cassette before covering its cover.
2. Plug one finger into rubber hole of the Oximeter (it is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.
3. Press button on the front panel;
【Note: if equipped with automatic startup function refers to clamp oximeter, need not press the button, the instrument has 5 s automatic signal detection function, directly inserted into the finger, instrument is automatically switched on timely】
4. Don't tremble your finger when the Oximeter is working. Your body is not recommended on moving status.
5. Press the button on the front panel, if we want change display direction;
【Note: if equipped with gravity sensing function of instrument then don't press the button, hand movements, the instrument with the gravity sensing has four corresponding interface switch】
6. Read relevant datum from display screen.
7. The instrument has the function of sleep, no signal 8seconds will enter standby state of sleep;
8. Please replace new batteries when OLED indicates the batteries are in low power.



Note: When plugging your finger into the Oximeter, your nail surface must be upward.

Declaration: Please use the medical alcohol to clean the rubber before each test and clean the tested finger with alcohol before and after the test. (The rubber inside of the Oximeter adopts medical rubber, which has no toxin, no harm, and brings no side effect such as allergy to the our skin).

Battery installation


1.According to the positive and negative identity right will be two AAA batteries into the battery door.

2.Push the battery cover horizontally in the direction of the arrow at the bottom.

Please note: pay attention to the battery positive and negative polarity, must be installed correctly, otherwise it may cause damage to instrument.


Lanyard installation

1, Put the lanyard thin end through the hole.



2, Put the lanyard coarser end through its already wearing thin end part and tighten.

Brief Description of Front Panel



Keys function description: in standby mode, start the key instrument into the working state, push down this button under working state, can change the display mode.

⚠Note⚠: The instruments equipped with gravity sensing function , OLED display mode only1. 2.3.4, figure 5 and 6 are not shown

⚠Note⚠: The machine profile picture only for your reference use, specific in kind prevail

Parameter setup:


Press start button (>3s) , into parameter setup .As menu ①:

1. In menu ①, When the “*” signal is shown on the “Alm Setup”, press the button (>3s) and enter into the menu ②. Press the button (<1s) can select item, then press button (>3s) to set the on/off for prompt ,beep,demo and screen brightness adjustment (optional "1" ," 2 "," 3 "and" 4 "). When the “*” signal is shown on the “Restore”,press the button (>3s) and all the settings are back to the factory settings.

⚠Attention⚠


- Using devices outdoors or under strong light, please adjust the screen brightness to a higher proper level for observation.
- It's better for user to choose a lower brightness to conserve battery power.

2. In menu ②, When the “*”signal is shown on the “Sounds Setup”, press the button (>3s) and enter into the menu ②,Press the button (<1s) can select item ,then press the button (>3s) to setup data.choose “+”or“-” to plus or minus values.





Menu ①



Menu ②

Note: If the matching refers to clamp oximeter is basic, not equipped with parameter setting function

Detailed descriptions of product functions:

1. Display Type: OLED display
2. SpO2: Measurement range: 70%~100%
 Accuracy: 80%~100%:±2% (Including 80%);
 70%~79%: ±3% ;
 Below 70% no requirement;
 Resolution: 1%
3. PR: Measurement range: 30BPM~254BPM
 Accuracy: ≤100BPM ,±1BPM
 >100BPM ,±2BPM
4. Parameters of LED probe

	wavelength	radiation power
RED	660±5nm	1.8mW
IR	940±10nm	2.0mW

This parameter can be especially useful to clinicians.


5. PI index: measurement scope: 0~20 (optional)
6. Power: two AAA 1.5V batteries
7. Automatic standby: the product shuts off by itself when no finger is in the product about 8 seconds
8. Automatic start up: every 5 s instrument will automatically detect the signal, after the hole with my finger, timely automatically boot;(optional)
9. Gravity sensing function: finger movement, the screen display will change with the gravity sensing changes (optional)
10. Dimension: about 58mm×36mm×33mm

11. Operation environment:

Temperature: 5 °C ~ 40 °C
 Humidity: 15% ~ 80%,no condensing
 Atmospheric pressure: 70kPa ~ 106kPa

Transport, storage environment:

Temperature: -10 °C ~ 40 °C
 Humidity: 10%~ 95%,no condensing
 Atmospheric pressure: 70kPa ~ 106kPa, non-corrosive gas and well-ventilated




environment.

12. Declaration: EMC of this product comply with IEC60601-1-2 standard.
13. This equipment has no alarm function.
14. Data averaging and signal processing delay the display and transmission of SpO2 data.
The date update period :<12S
15. About 30min required for oximeter to warm from the minimum storage temperature between uses until it is ready for intended use
16. About 30min required for oximeter to cool from the maximum storage temperature between uses until it is ready for intended use
17. Applied parts specified: Probe and its circuit
18. Use specification

Expected medical instructions	The fingertip pulse oximeter can be used to measure human haemoglobin saturation and pulse rate through finger
Expected patient population	-People who need blood oxygen measurement .
Expected use or interaction with body parts tissue type	Finger
Expected user profile	People who need blood oxygen measurement , doctors, etc
Application environment	Avoid electromagnetic interference Extreme temperature Avoid pollution and dust Avoid direct sunlight, etc
Operating principle	Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 940nm near infrared light) can be focused onto human nail through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on two groups of diodes through process in electronic circuits and microprocessor.

Classification

1. Management Class for Medical Devices: IIequipment
2. Anti-electric Shock Type: Internally powered equipment
3. Anti-electric Shock Degree: Type BF equipment
4. Overvoltage category classification:Class I
5. Pollution degree:



Pollution degree2:Micro-environment with non-conductive pollution,expect occasional conductivity caused by condensation

Maintenance and Preservation

1. Replace the batteries timely when the low battery indicator flashes
2. Clean the surface of fingertip pulse oximeter before use.
3. Remove the batteries when the oximeter is not likely to be used for some time when leakage from batteries would result in an unacceptable risk.
4. It would be better to preserve the product in -10~40°C (14-104°F) and humidity is 10%-95%.
5. It is recommended that the product should be kept dry anytime. A wet ambience might affect its lifetime and even damage the product.

6. Cleaning frequency

If the oximeter is dirty when used at home, it is recommended to clean the enclosure and silicone pad after each use.If the oximeter is not dirty, simply clean the black silicone pad before and after each use.

When used in a medical institutions, clean it after each use.

7. Cleaning method

When the oximeter is used at home, wipe and disinfect with 75% alcohol, then dry naturally or clean the product with dry cloth. Clean at least twice a week.

When using in medical institutions, clean the tested finger with 75% alcohol before use. After each measurement, wipe and disinfect the contact part between black silicone and enclosure with 75% alcohol.

8. The oximeter can display functional arterial oxygen saturation and pulse rate after calibration.
9. The name of the simulator is FLUKE Index2 simulator, and the version number is 3.00.
10. The Blood oxygen simulator is calibrated to display functional oxygen saturation.
11. The maximum temperature of the contact surface between the product and human body does not exceed 41 °C.

Waste disposal

- 1.Please follow local laws to dispose of waste scrap.
2. Follow local ordinances and recycling instructions regarding to disposal or recycling of the device and device components, including used batteries and packaging box


Product Accessories

Inventory list	Quantity
The host	1 set
Lanyard	1 pcs
One user manual	1 pcs
AAA battery (optional)	2pcs

Manufacturer's Declaration of the EUT

Statement:

1. The Fingertip Pulse Oximeter or user should use the product in the electromagnetic environment specified in the following table, otherwise it may cause abnormal operation of the



product.

2. Fingertip Pulse Oximeter is a table-top equipment, it suitable for medical unit and home use.
3. Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm(12 inches) to any part of the oximeter, including cables specified by the manufacturer. Otherwise, degradation of the performance of this Fingertip Pulse Oximeter could result
4. Warning: Use of this Fingertip Pulse Oximeter adjacent to or stacked with other equipment should be avoided because it could result in improper observed to verify that they are operating normally.
5. Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.


Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

1	Guidance and manufacturer's declaration – electromagnetic emission	
2	The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of Fingertip Pulse Oximeter should assure that it is used in such an environment.	
3	Emissions test	Compliance
4	RF emissions CISPR 11	Group 1
5	RF emissions CISPR 11	Class B
6	Harmonic emissions IEC 61000-3-2	N/A
7	Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A


Guidance and manufacturer's declaration – electromagnetic immunity –

for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity
The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Pulse Oximeter should assure that it is used in



such an environment.		
Immunity test	EN 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15kV air
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	N/A
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines	< 5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T)	N/A




IEC 61000-4-11	<p>for 5 cycles</p> <p>70 % U_T (30 % dip in U_T)</p> <p>for 25 cycles</p> <p>< 5 % U_T (>95 % dip in U_T)</p> <p>for 5 sec</p>	
<p>Power frequency</p> <p>(50/60 Hz)</p> <p>magnetic field</p> <p>IEC 61000-4-8</p>	30 A/m	30A/m

Guidance and manufacturer's declaration – electromagnetic immunity –

for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity		
The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Pulse Oximeter should assure that it is used in such an environment.		
Immunity test	EN 60601 test level	Compliance level



Conducted RF	3 Vrms	N/A
IEC 61000-4-6	150 kHz to 80 MHz	

Table 9 - Test specifications for enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL
385	380 -390	TETRA 400	Pulse modulation ^{b)}	1,8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0.3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
870						
930						
1 720	1 700-1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1 845						
1 970						

Xuzhou Yongkang Electronic Science Technology Co., Ltd.
 Only for _____ Registration _____ to use
 12

2 450	2 400-9 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5 240	5 100-5 800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0.3	9
5 500						
5 785						

NOTE:

If necessary to achieve the immunity test level, the distance between the transmitting antenna and the me equipment or me system may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.


c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.


















Possible Problems and Resolutions

Problem	Possible reason	Solution
SpO2 or PR can not be shown normally	1. Finger is not plugged correctly 2. Patient's Oxyhemoglobin value is too low to be measured	1. Retry by plugging the finger 2. Try more times. If you can make sure there is no problem in the product, please go to hospital timely for exact diagnosis
SpO2 or PR is shown unsteady	1. The finger might not be plugged deep enough 2. Finger is trembling or the patient is on movement status	1. Retry by plugging the finger 2. Please remain at rest
The Oximeter can not be turned on	1. Inadequate power or power off 2. Batteries might be installed incorrectly 3. The Oximeter might be damaged	1. Please replace the batteries 2. Please reinstall the batteries 3. Please contact with local customer service centre
Indication lamps are suddenly off	1. The product automatically shuts off when no signal is detected in 8 seconds 2. Inadequate power	1. Normal 2. Replace the batteries

Symbols and Definitions

	BF type application part		Serial number
---	--------------------------	---	---------------




	To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.		Date of manufacture
	Reference manual		Manufacturer
	Cautions		European union representative
IP22	First characteristic numeral 2: Against ingress of solid foreign objects:≥12.5mm diameter Second characteristic numeral 2 Against ingress of water with harmful effects:dripping(15° tilted)		Standby
	Humidity range		Temperature range
	Keep dry		Avoid sunlight
	Product certification		Up toward
	Lot number		This device has no alarm function
	signal inadequacy		①Indication of probe faults (open circuit condition or close circuit condition) ②Indication of Probe cable faults ③Indication of Probe cable extender faults

Reserves the right to technical change appearance, our products are subject to change without prior notice, please forgive me!

Statement:

1, If you need maintenance, please contact the manufacturer

2, The company can be in the form of email or other electronic files provide users with random files.



3, The instrument is not used for evaluation of blood oxygen probe pulse and pulse blood oxygen monitor accuracy.

After-sales service

Ensure that users

- Please read user manual before using the instrument;
- According to the requirement of the instruction manual for the operation and daily maintenance, and make sure the machine power supply, and environmental requirements

Maintenance regulations

- To conform to the regulations, free maintenance within the scope of products, with warranty card for free maintenance. All that is beyond the scope of free maintenance product, provide paid services.
- With warranty card and shopping invoice, main machine for a year, accessories for three months are under free maintenance services from the date of purchase.
- Following does not belong to the scope of free maintenance
 - ☞ The fault caused by human factors, the damage;
 - ☞ Due to the use to be inconsistent with the provisions of our company work environment to cause damage;
 - ☞ Due to the product in the our company authorized personnel disassembling or repairing damaged;
 - ☞ Products beyond the warranty period.

Maintenance Time

- If any problem, please call us in 9: 00 am to 5: 30 pm from Monday to Friday(except national holiday), call us: 400-828-6667

Company information


Using life: 5 years

TEL: +86-516-87892766 -601 FAX: +86-516-87892755-606


Web site: www.yonker.cn

Date: see product label

Xuzhou Yongkang Electronic Science Technology Co., Ltd

 1st&2nd Floor, 6#01, 6#02, No.6 Building 1st Phase Economic Development Manufacturing Zone, LIANDO U Valley, No.6 Leye Road, Xuzhou ETDZ, 221000 Xuzhou, PEOPLE'S REPUBLIC OF CHINA

The authorized representative of European Union:



Prolinx GmbH
Brehmstr. 56, 40239 Düsseldorf
Germany

CE 0123

