5. Tip-to-Tip Pinch

Place the tip of the thumb on the center of the ergonomic palm bar. Position the finger to be exercised on one of the center buttons. Keeping both tips slightly bent, flex the thumb and the finger being exercised toward the center of the unit at the same time.

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6. Fingertip Flexion

Place all finger tips on their own button, with the tip of the thumb on the center of the ergonomic base. Keeping all finger tips bent, flex the thumb and all the fingers toward the center of the unit at the same time.

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
PROGRESSIVE EXERCISE SYSTEMS

DIGI-FLEX®

exercise each finger individually, or entire hand

Resistance level: per finger / per hand
Weerstandsniveau: per vinger / geheel
Niveau de résistance: par doigt / sur l'ensemble
Widerstandskraft: pro finger / gesamtes gerät
Livello di resistenza: per singolo dito / intero apparecchio
Nivel de resistencia: por dedo / sistema completa

10-0751	10-0740	10-0741	10-0742
tan xx-light braun xx-leicht beige xx-licht beige xx-légère marrone xx-leggero marfil xx-ligero	yellow x-light gelb x-leicht geel x-licht jaune x-légère giallo x-leggero amarillo x-ligero	red light rot leicht rood licht rouge légère rosso leggero rojo ligero	green medium grün mittel groen medium vert moyen verde medio verde mediano
0.75 lbs / 2.5 lbs 0.34 kgs / 1.1 kgs	1.5 lbs / 5.0 lbs 0.7 kgs / 2.3 kgs	3.0 lbs / 10.0 lbs 1.4 kgs / 4.5 kgs	5.0 lbs / 16.0 lbs 2.3 kgs / 7.3 kgs




10-0743	10-0744	10-0752	10-0753
blue heavy blau stark blauw zwaar bleu fort blu pesante azul pesado	black x-heavy schwarz x-stark zwart x-zwaar noir x-fort nero x-pesante negro x-pesado	silver xx-heavy silber xx-stark zilver xx-zwaar argent xx-fort argento xx-pesante plata xx-pesado	gold xxx-heavy golden xxx-stark goud xxx-zwaar doré xxx-fort oro xxx-pesante oro xxx-pesado
7.0 lbs / 23.0 lbs 3.2 kgs / 10.4 kgs	9.0 lbs / 30.0 lbs 4.1 kgs / 13.6 kgs	11.0 lbs / 36.0 lbs 5.0 kgs / 16.3 kgs	13.0 lbs / 43.0 lbs 5.9 kgs / 19.5 kgs

now available:

Cando® Digi-Extend® exerciser for finger extension (a)

Cando® Digi-Squeeze® finger exerciser for hand grasp and finger flexion (b)

Cando® Digi-Extend n'Squeeze® for flexion and extension combination exercises (c)

FEI
FABRICATION
ENTERPRISES INC

Fabrication Enterprises, Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523

tel: 800-431-2830 • +1-914-345-9300
fax: 800-634-5370 • +1-914-345-9800
FabEnt.com

Cando®, Digi-Flex® are trademarks of Goldberg
©Copyright FEI, 2020 - all rights reserved

CE

Authorized CE Representative:
AJW Technology Consulting GmbH
Königsallee 106
40215 Düsseldorf (Germany)

Digi-Flex® Exercise System:

10-0751, 10-0740, 10-0741, 10-0742, 10-0743, 10-0744, 10-0752, 10-0753

CE

INDICATIONS

- Develops individual finger strength, flexibility and coordination
- Also builds hand and forearm strength

PRECAUTIONS

- **If the individual is being treated for a hand disorder, specific exercises and the intensity and number of repetitions should be approved by a healthcare professional.**
- **Following injury or surgery, consult a healthcare professional before performing exercises.**

- **If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**

GUIDELINES FOR USE

Exercises should be done at a slow, comfortable pace. Suggested exercises follow.

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1. Gross Grasp

Place the middle of the fingers on the buttons with the finger tips wrapped over the top. Positioning the hook of the ergonomic palm bar over the web space between the thumb and index finger, flex the thumb and all fingers toward the center of the unit as if making a fist.

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
2. Hook Position

Place finger tips on individual buttons and position ergonomic palm bar with the hook over the web space between the thumb and index finger. Keeping finger tips slightly bent, flex all fingers toward the center of the unit.

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
3. Key Pinch

Place the hook of the ergonomic palm bar on the inside of the bent index finger. With the tip of the thumb slightly bent, flex the thumb on the first button.

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
4. Trigger Pinch

Loosely hold the unit vertically with the buttons toward the fingers and the hook of the ergonomic palm bar positioned over the web space between the thumb and index finger. With the middle of the index finger, press the second button as if squeezing a trigger.

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