


BASELINE[®] EVALUATION INSTRUMENTS

plastic 30 cm 360° HiRes[®] goniometer

REF

12-1000HR



- permits observation of joint's axis of motion, and range of motion
- 360° head has three scales for ease of use
- primary scale reads in 1° increments
- goniometer arm has a linear scale in inches and centimeters
- goniometer arm remains at measurement until reset

other goniometer styles and sizes available:

30 cm 360°	12-1000*
20 cm 360°	12-1001*
15 cm 360°	12-1002*
15 cm 180° pocket.....	12-1005*
15 cm 360° rulongmeter	12-1006*
7.6 cm -30° to 120° finger.....	12-1012*

* available standard or with white HiRes[®] strip behind numbers (add HR to end of product code for HiRes[®] model)



FEI
FABRICATION
ENTERPRISES INC.

Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300
fax: +1-914-345-9800
FabEnt.com

MD

EC REP



Authorized CE representative:
AJW Technology Consulting GmbH
Breite Strasse 3
40213 Dusseldorf Germany

ver 0122