

# **BASELINE<sup>®</sup>** EVALUATION INSTRUMENTS

plastic 20 cm 360° HiRes<sup>®</sup> goniometer

**REF**

**12-1001HR**



- permits observation of joint's axis of motion, and range of motion
- 360° head has three scales for ease of use
- primary scale reads in 1° increments
- goniometer arm has a linear scale in inches and centimeters
- goniometer arm remains at measurement until reset

## **other goniometer styles and sizes available:**

30 cm 360° .....	12-1000*
20 cm 360° .....	12-1001*
15 cm 360° .....	12-1002*
15 cm 180° pocket.....	12-1005*
15 cm 360° rulongmeter .....	12-1006*
7.6 cm -30° to 120° finger.....	12-1012*

\* available standard or with white HiRes<sup>®</sup> strip behind numbers (add HR to end of product code for HiRes<sup>®</sup> model)



**FEI**  
FABRICATION  
ENTERPRISES INC.

**EC REP**



Fabrication Enterprises Inc.  
250 Clearbrook Rd, Suite 240  
Elmsford, NY 10523 (USA)  
tel: +1-914-345-9300  
fax: +1-914-345-9800  
FabEnt.com



Authorized CE representative:  
AJW Technology Consulting GmbH  
Breite Strasse 3  
40213 Dusseldorf Germany

ver 0122