


# **BASELINE<sup>®</sup>** **EVALUATION INSTRUMENTS**

**REF** 12-0297

## **PNEUMATIC BULB HAND DYNAMOMETER / PINCH GAUGE**



Adjustable dynamometer/pinch gauge measures both grip and pinch strength. Comes with three different sized interchangeable squeeze bulbs that can be easily fitted to the gauge. Lightweight and portable. Maximum reading remains until reset by turning center knob. Unit is calibrated in PSI (0-30) and kPa (0-200).

**⚠ This Product Contains Dry Natural Rubber.**

**FEI** Fabrication Enterprises Inc.  
250 Clearbrook Rd, Suite 240  
Elmsford, NY 10523 (USA)  
tel: +1-914-345-9300 • 800-431-2830  
fax: +1-914-345-9800 • 800-634-5370  
FabEnt.com

**EC REP**  
**CE** AJW Technology Consulting GmbH  
Breite Strasse 3  
40213 Düsseldorf (Germany)

**MD**