


BASELINE[®] EVALUATION INSTRUMENTS

Pneumatic Bulb Dynamometer with reset



12-0291

5 cm squeeze bulb dynamometer measures PSI.

Red maximum force indicator holds maximum force reading until reset. Max indicator can be reset to zero by turning the center reset knob counterclockwise.


12-0291: 0-30 PSI, with max pointer, 0.5 PSI gradations

12-0293: 0-15 PSI, with max pointer, 0.5 PSI gradations

This product contains dry natural rubber.



Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com



AJW Technology Consulting GmbH
Breite Strasse 3
40213 Düsseldorf (Germany)

